DISINFECTANT CLEANER & DEODORIZER
WIPES

Reference Sheet

Manufactured by:
Virox Technologies Inc.
2770 Coventry Road
Oakville, ON Canada L6H 6R1
1-800-387-7578 | virox.com
viroxanimalhealth.com

FOR VETERINARY USE

For use in animal health facilities and areas. For disinfection and cleaning of hard, non-porous surfaces found in animal premises facilities such as veterinary clinics and offices, animal shelters, kennels, treatment facilities, transportation vehicles, pet grooming facilities and pet shops.

VIRUCIDAL IN 1 MINUTE

Adenovirus Type 5
Adenovirus Type 8
Avian Influenza A Virus (H3N2)
Avian Influenza A Virus (H7N9)
Avian Reovirus
Canine Distemper Virus
Canine Influenza A Virus (H3N2)
Canine Parvovirus (CPV) (Strain Cornell)
Cytomegalovirus
Feline Calicivirus
Feline Herpes Virus (Feline Rhinotracheitis)
Feline Infectious Peritonitis Virus (FIPV)
Feline Leukemia Virus (FeLV)
Feline Panleukopenia Virus
Hepatitis B Virus (HBV)
Hepatitis C Virus (HCV)
Herpes Simplex Virus Type 1
Herpes Simplex Virus Type 2
HIV Type 1 (Strain HTLV III_B)
Human Coronavirus
Influenza A Virus (Strain A/8/68)
Influenza B Virus
Newcastle Disease Virus
Norovirus (Feline Calicivirus as the surrogate)
Pandemic Influenza A Virus H1N1
Parainfluenza Virus Type 3
Poliovirus Type 1 (Strain Chat)
Pseudorabies Virus
Respiratory Syncytial Virus (RSV)
Rhinovirus Type-37
Rotavirus (WA strain)
SARS-CoV-2 (the virus that causes COVID-19)
(SARS-Related Coronavirus 2)
Vaccinia Virus

This product has demonstrated effectiveness against SARS-CoV-2 (the virus that causes COVID-19) (SARS-Related Coronavirus 2).

Proven effectiveness against Adenovirus Type 5 which allows for a Broad-Spectrum Virucide claim against most enveloped and non-enveloped viruses.

BACTERICIDAL IN 1 MINUTE

Acinetobacter baumannii (Multi-Drug Resistant (MDR))
Bordetella bronchiseptica
Bordetella pertussis
Burkholderia cepacia
Enterobacter aerogenes
Enterobacter cloacae (New Delhi metallo- β -lactamase 1 (NDM -1))
Enterococcus faecalis (Vancomycin Resistant (VRE))
Enterococcus faecium
Escherichia coli
Escherichia coli O157:H7
Escherichia coli (Extended Spectrum β -lactamase (ESBL))
Escherichia coli (New Delhi metallo- β -lactamase 1 (NDM-1))
Klebsiella pneumoniae
Klebsiella pneumoniae (Carbapenem Resistant (CRE))
Klebsiella pneumoniae (Extended Spectrum β -lactamase (ESBL))
Klebsiella pneumoniae (New Delhi metallo- β -lactamase 1 (NDM-1))
Legionella pneumophila
Listeria monocytogenes
Pasteurella multocida
Pseudomonas aeruginosa
Pseudomonas aeruginosa (Metallo- β -lactamase (MBL) positive)
Salmonella enterica (formerly known as *Salmonella choleraesuis*)
Serratia marcescens
Shigella dysenteriae
Staphylococcus aureus
Staphylococcus aureus (Community Acquired Methicillin Resistant Genotype USA 400)
Staphylococcus aureus (Methicillin Resistant (MRSA))
Staphylococcus aureus (Vancomycin Intermediate Resistant (VISA))
Staphylococcus aureus (Vancomycin Resistant (VRSA))
Staphylococcus epidermidis (Methicillin Resistant (MRSE))
Staphylococcus pseudintermedius
Stenotrophomonas (Pseudomonas) maltophilia
Streptococcus pneumoniae (Penicillin Resistant (PRSP))
Streptococcus pyogenes

FUNGICIDAL IN 1 MINUTE

Candida albicans
Microsporum canis
Trichophyton interdigitale (formerly known as *Trichophyton mentagrophytes*)

MILDEWICIDAL IN 4 MINUTES

Aspergillus brasiliensis

TUBERCULOCIDAL IN 1 MINUTE

Mycobacterium bovis (BCG)

BROAD-SPECTRUM NON-FOOD CONTACT SANITIZER IN 10 SECONDS

Acinetobacter baumannii

Campylobacter jejuni

Enterobacter aerogenes

Enterococcus faecalis (Vancomycin Resistant (VRE))

Escherichia coli O157:H7

Listeria monocytogenes

Pseudomonas aeruginosa

Salmonella enterica (formerly known as *Salmonella choleraesuis*)

Staphylococcus aureus

Staphylococcus aureus (Methicillin Resistant (MRSA))

DISINFECTANT USE DIRECTIONS

HEAVILY SOILED SURFACES REQUIRE CLEANING PRIOR TO DISINFECTION

For Use as a One-Step Cleaner/Disinfectant:

1. Pre-clean visibly soiled areas.
2. Pull towelette from dispenser and wipe hard, non-porous surfaces.
3. Allow surfaces to remain visibly wet for 1 minute.
4. Allow to air dry.

To Disinfect Non-Critical, Pre-Cleaned Instruments:

1. Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and rough dried.
2. Wipe all hard non-porous surfaces of instruments with this product until thoroughly wet.
3. Allow surfaces to remain visibly wet for 1 minute.
4. Wipe surfaces dry. If streaking is observed, wipe with a clean, damp cloth.

Animal Housing Facilities in Veterinary Clinics and Hospitals:

For cleaning and disinfecting the following hard, non-porous surfaces: feeding and watering equipment, instruments, cages, kennels, stables, catteries, etc.

1. Remove all animals and feeds from premises, animal transportation vehicles, crates, etc.
2. Remove all litter from floors, walls and surfaces of facilities occupied by animals.
3. Empty all feeding or watering appliances.
4. Thoroughly clean all surfaces with soap and detergent and rinse with potable water. Pull towelette from canister and wipe hard, non-porous environmental surfaces.
5. All surfaces must remain visibly wet for 1 minute.
6. Allow to air dry.*
7. Do not house animals or employ equipment until treatment has been absorbed, set or dried.
8. Thoroughly scrub all treated automatic feeders and waterers with soap or detergent, and rinse with potable water before reuse.

To Kill Mold and Mildew:

1. Pre-clean surfaces.
2. Pull towelette from dispenser and wipe product onto hard, non-porous surfaces until completely wet.
3. Allow the surfaces to remain visibly wet for 4 minutes.
4. Allow to air dry.

To Control Mold and Mildew:

1. Pull towelette from dispenser and wipe product onto hard, non-porous surfaces until completely wet.
2. Allow to air dry.
3. Repeat application weekly or when growth reappears.

Broad-Spectrum Sanitizing on Environmental Surfaces:

1. Pre-clean visibly soiled areas.
2. Pull towelette from dispenser and wipe hard, non-porous surfaces.
3. Allow surface to remain wet for 10 seconds.
4. Wipe dry. No rinse required. Food contact surfaces require rinsing with potable water.

For Use as a Cleaner and/or Deodorizer:

1. Pull towelette from dispenser and wipe hard, non-porous environmental surfaces.
2. Allow to air dry.

PRECAUTIONARY STATEMENTS

Avoid (i) contact with eyes; (ii) contamination of food in the application and storage of this product and; (iii) mixing with other cleaning or disinfecting products.

FIRST AID

If in contact with eyes, flush immediately and thoroughly with water. Call a physician if irritation persists. If ingested, call a physician or a poison control center immediately and do not induce vomiting. TAKE LABEL OR PRODUCT NAME AND DIN WITH YOU WHEN SEEKING MEDICAL ATTENTION.

STORAGE

Avoid storage at elevated temperatures.

DISPOSAL

For information on the disposal of unused, unwanted product and the clean-up of spills, contact the Provincial Regulatory Agency or the Manufacturer. Rinse the container with water. Dispose of container in accordance with municipal, provincial and federal regulations.

REFILLABLE CONTAINER

Refill this container with this brand name only.

DISCLAIMER: You must read and follow the instructions on the label for our product and use the product for the particular applications specified on the product label. To the full extent permitted under applicable law, Virox Technologies Inc. shall not be liable for any damages whatsoever, arising out of or otherwise in connection with any failure to follow those instructions or any use of its products which is not specified on the product's label. The instructions on the product label are not a warranty, express or implied, and Virox Technologies, Inc. expressly disclaims any implied warranties, including any warranty of fitness for a particular purpose. Under applicable law, some of the above waivers or disclaimers of warranty, or exclusions of liability for certain types of damages may not apply to you. In such cases, our liability will be limited to the fullest extent permitted by applicable law.

For relevant patents, visit: www.virox.com/patents

Accelerated Hydrogen Peroxide® and Design are registered trademarks of Diversey Inc.

* Ventilate buildings and other closed spaces.

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Accelerated Hydrogen Peroxide® et le dessin sont des marques déposées de Diversey Inc.